

DETROSAN SFC – CONCENTRATE SURFACE DISINFECTANT

DESCRIPTION OF PRODUCT:

QAC and Amine based floor and surface disinfectant

SPECIFICATION OF PRODUCTS:

Chemical-Physical Characteristic

Parameters	M.U	Specifications
Physical form	-	Liquid
Appearance	-	Limpid
Colour	-	Clear
Odour	-	Characteristic
pH	pH unit, 20'C	9-11
Density (20°C)	g/cm ³ , 20'C	1,04 ± 0,1 g/cm ³
Solubility in water	-	Soluble
Benzalconium Chloride	W/W	7,5%±%15
DDAC	W/W	5%±%15
Dodecylamine	W/W	2%±%15

*Confirmed by Detrox Quality Control Department.

Microbiological Characteristics

Parameter	M.U.	Specifications	Methods of Analysis
Escherichia coli NCTC 10538	cfu/ml	%0,25 concentration, Effective in 5 min.	EN 1276
Staphylococcus aureus ATCC 6538	cfu/ml	%0,25 concentration, Effective in 5 min.	EN 1276
Pseudomonas aeruginosa ATCC 15442	cfu/ml	%0,25 concentration, Effective in 5 min.	EN 1276
Enterococcus hirae ATCC 10541	cfu/ml	%0,25 concentration, Effective in 5 min.	EN 1276
Mycobacterium avium ATCC 15769	cfu/ml	%0,25 concentration, Effective in 15 min.	EN 14348
Mycobacterium terrae ATCC 15755(Strain, W45)	cfu/ml	%0,25 concentration, Effective in 15 min.	EN 14348
Candida albicans ATCC 10231	cfu/ml	%0,25 concentration, Effective in 15 min.	EN 1650

DETROSAN SFC – CONCENTRATE SURFACE DISINFECTANT

Aspergillus niger ATCC16404	cfu/ml	%0,25 concentration, Effective in 15 min.	EN 1650
Adeno virüs ATCC VR-5	cfu/ml	%0,25 concentration, Effective in 15 min.	EN 14476
Poliovirüs ATCC VR-192	cfu/ml	%0,25 concentration, Effective in 15 min.	EN 14476
Murine Norovirüs PTA-5935	cfu/ml	%0,25 concentration, Effective in 15 min.	EN 14476

INTENDED USE OF PRODUCT:

Usage Area

- It is used for cleaning and disinfection of surfaces existing at hospitals and health foundations, like surgery rooms, clean rooms, intensive care units, patient rooms, new born units etc.

Ingredients

- Bis (3-aminopropyl) dodecylamine, Didecyldimethylammonium chloride, Benzalkonium chloride

Usage form

- Detrosan SFC is diluted at %0,25 concentration (2,5 ml product into 1 L water)
- Floor or medical device surfaces are cleaned with moistened wipe or mop with Detrosan SFC.
- Wait for 15 min. effect time.
- Life of prepared solution is 24 hours. In the case of intensive contamination, it is recommended to change solution.

Microbiological effect

- Bactericidal (incl.Tbc), fungicidal, virucidal effect in 15 minutes - %0,25 concentration

First aid Informations

- Eyes; In the case of eye contact, wash plenty of water around 5 min.
- Skin: Wash your skin immediately.
- Ingestion: Drink water and apply your doctor.

Netice TONBİL
General Manager

DETRO HEALTHCARE KİMYA SAN. A.Ş.
Mavi Sokak Kat:1 Kat:1 Kat:1 Kat:1 Kat:1
No:1 Kat:1 Kat:1 Kat:1 Kat:1 Kat:1
Tel:0212 659 77 62 - 0212 659 77 63
Fax:0212 659 77 63 - 0212 659 77 63