

# BACTINYL LINGETTES DESINFECTANTES PARFUMÉES SINGLE USE PERFUMED DISINFECTANT WIPES



## BENEFITS

- Safe for surfaces even the most fragile
- Bactericidal, yeasticidal and virucidal
- Food contact agreement



Formule douce



## APPLICATION

Single use wipes impregnated with disinfectant solution for cleaning and disinfection by contact of non-invasive medical devices, instrumentation, medical equipment and equipment.

Medical device. CLASS IIA.

## INSTRUCTIONS FOR USE

Rub the surfaces for with BACTINYL LINGETTES DESINFECTANTES PARFUMÉES. Respect the action time, after the passage of the wipe, according to the efficiency desired. Let dry. Do not rinse except in case of contact with surfaces that may come into contact with foodstuffs.

Close the box after use for the impregnation rate is optimum to the last wipe. Do not replace disinfection by soaking.

Complies with the decree of 19/12/13 amending the decree of 08/09/1999 concerning products for cleaning materials and articles intended to come into contact with foodstuffs.

## CHARACTERISTICS OF THE IMPREGNATING LIQUID

Appearance :	liquid
Colour :	colourless
Perfume :	pine
Density :	= 1
pH :	low base
Expiration :	2 years
Stockage :	4°C - 35°C



## MICROBIOLOGICAL ACTIVITY OF THE IMPREGNATING LIQUID

**Bactericidal :** EN 13727 (dirtiness conditions – 5 min – 20°C), EN 14561 (dirtiness conditions – 5 min – 20°C), EN 16615 (dirtiness conditions – 2 min – 20°C)

**Yeasticidal :** EN 13624 (dirtiness conditions – 5 min – 20°C), EN 14562 (dirtiness conditions – 15 min – 20°C), EN 16615 (dirtiness conditions – 2 min – 20°C)

**Sporicidal :** EN 13704 (cleanliness conditions – 60 min – 20°C)

**Virucidal :** EN 14476 on PRV\*, BVDV\*, H1N1, Herpès, Vaccinia virus, Rotavirus (dirtiness conditions – 5 min – 20°C), Norovirus (dirtiness conditions – 15 min – 20°C)

\* Swine Pseudorabies Virus (PRV) and Bovine Viral Diarrhea Virus (BVDV) are surrogate of HBV and HCV.

## PACKAGING

Tub of 120 wipes

Réf. LB920600

Tub of 120 wipes (large size)

Réf. LB920740

## SECURITY

For professional use only.

For more information please refer to the MSDS.

